




Ammonia Chlorine Solution for EZ1102 Ammonia Analyser, 500 mL

Kat. číslo produktu:	APPC1102-01
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Ready-to-use Reagents for EZ Series Online Analysers

Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

Technické údaje

Obsah balení:	500 mL
Parametr:	Ammonia