



## Test kit, nitrogen, nitrate color disc, model NI-14

**Kat. číslo produktu:**  
CZK Cena bez DPH:  
Datum odeslání neuvedeno

**1416100-EU**  
Kontaktujte nás

Nebezpečné



Utilizing the simple color disc method, this kit measures Nitrogen, Nitrate as NO<sub>3</sub> and contains 50 tests.

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Technické údaje

Metoda:	Color disc/ Cadmium reduction
Nejmenší dílky přírůstku:	0,05
Parameter:	Nitrogen, Nitrate
Parametr:	Nitrogen, Nitrate - As NO <sub>3</sub> <sup>-</sup> -N
Platforma:	Barevný disk
Počet testů:	50
Přístroj:	NI-14
Rozsah měření:	Dual Range
Rozsah měření (2):	0 - 10 mg/L NO <sub>3</sub> <sup>-</sup> -N
Smallest increment:	0,05
	0,5
Standardní rozsah měření:	0 - 1 mg/L NO <sub>3</sub> <sup>-</sup> -N; 0 - 10 mg/L NO <sub>3</sub> <sup>-</sup> -N
Styl kufříku:	B
Typ testovací soupravy:	NI-14**