



Be Right™



## Nitrate test kit, model NI-11

Kat. číslo produktu:

146803-EU

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno

Nebezpečné



One of Hach's most popular test kits, it uses the simple, low-cost color disc method.

**One of Hach's most popular test kits**

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Technické údaje

Metoda:	Color disc/ Cadmium reduction
Parameter:	Nitrogen, Nitrate
Parametr:	Nitrogen, Nitrate - As $\text{NO}_3^-$ -N
Platforma:	Barevný disk
Počet testů:	100
Přístroj:	NI-11
Rozsah měření:	Single Range
Smallest increment:	1.0
Standardní rozsah měření:	0 - 40 mg/L $\text{NO}_3^-$ -N
Styl kufříku:	B
Typ testovací soupravy:	NI-11**